

15 mL Digestion Vessels and Accessories



Our 15 mL digestion vessels (MP-123 and MP-123PN) are certified to class A volume accuracy, include a metals background report and are compatible with leading brand graphite block digestion systems.

- 1,000 vessels/case
- Made from virgin polypropylene homopolymer (PP-H)
- Excellent visibility
- Very good resistance to acids
- Molded-in graduations are in 1 mL increments to 15 mL
- 5, 10- and 15-mL graduation lines extend around the entire circumference
- Each shipment comes with a certificate of volume accuracy stating the tubes meet Class A tolerance at 5, 10- and 15-mL graduation marks
- Total capacity of 18 mL
- Extremely low leachable metal content; each shipment comes with a metals background report
- Fit leading brand graphite digestion block systems
- Can be used as a digestion vessel, graduated cylinder, volumetric flask, autosampler tube, or storage container
- Screw caps are available in natural
- We also offer ribbed watch glasses (MP-128)

Our MP-123 15 mL digestion tube is a direct replacement for Environmental Express SC415 and SCP Science 010-515-263 and fit the EE SC196 96 well HotBlock and 15ml DigiPREP Mini, Jr and LS systems.

MP-128 ribbed watch glass is equivalent to Environmental Express SC417.