

June 2, 2020

## 100 mL Titration Sample Beaker



MP-118 100 mL, 53.6 x 90.5mm is a widely used size for titration/TAN/TBN systems by Mettler-Toledo (M-T). Equivalent to M-T 00101974, 511093388 & 51109389; Dyn-a-med 80099; Precision Weighing Balances 99008; Fisher Scientific NC9485544; VWR 56221-531 & 12578-108.

- 450 beakers/case
- Made from virgin polypropylene
- Graduations at 20, 40, 60, 80 and 100 mL
- Excellent visibility
- Good resistance to acids

Fits the following M-T sample introduction systems:

- Rondolino™
- Rondo™
- InMotion™ Flex Autosampler
- InMotion™ Pro Autosampler
- InMotion™ Max Autosampler

Used on the following M-T instruments:

- G10S
- G20S
- T5
- T7
- T9