

May 20, 2020

120 mL Titration Sample Beaker



MP-117 120 mL, 43.3 x 113 mm is a common beaker size for various titration methods especially TAN and TBN. Equivalent to Metrohm 6.1459.300 and 6.1459.307 absent threads and threaded cap option. Fits multiple USB and Robotic Sample platforms used with Metrohm Titration systems.

- 300 beakers/case
- Made from virgin polypropylene
- Excellent visibility
- Good resistance to acids

Fit Metrohm sample introduction systems:

- 814 USB Sample Processor
- 815 Robotic USB Sample Processor XL
- 815 Robotic Soliprep
- 855 Robotic Titrosampler
- 864 Robotic Balance Sample Processor
- 898 XYZ Sample Changer
- MATi 02

Fit Metrohm racks:

- 6.2041.470
- 6.2041.840
- 6.2041.850
- 6.2068.010

May 20, 2020

- 6.2068.020
- 6.2071.020
- 6.2071.030
- 6.9920.141

Methods:

- ASTM D 664
- ASTM D 2896