

March 10, 2020

CAV Shell Vials



MoldPro automatic viscometer cups are cylindrical, translucent sample vials for Cannon CAV[®] Automatic Viscometers.

- Made from virgin polypropylene
- Easily withstands 15 – 100°C temperature range
- Used for ASTM D445/D446, ISO 3104/3105, IP 71
- Excellent visibility
- Good resistance to acids
- 1000 cups/case

March 10, 2020

Item	Volume	Dimension	Substitute or replacement for	Instruments	Sample Introduction
MP-104	15 mL	24 x 62 mm	Cannon® 61.3663 & 61.3661 (thin wall)	Cannon Automatic Viscometers (CAV®), CAV® 4.2 Dual-Bath Kinematic Viscometer, CAV®- 2100 & CAV®- 2200 Single – Bath Kinematic Viscometer	Cannon® 13 - & 14 - position sample handler