

March 26, 2020

Ion Chromatography (IC) Sample Vials



MoldPro sample vials are a direct replacement for Dionex items 055058 & 074228 and Environmental Express item K4300 for Dionex AS-50 and AS-AP autosamplers.



March 26, 2020

- Made from virgin polypropylene
- Pierce cap feature
- 500 vials and caps/case

Item	Volume	Dimension	Substitute or direct replacement for	Instruments	Sample Introduction/Autosampler	Trays/Racks
MP-113P	~7 mL	20 x 46 mm	Dionex™ 055058 & 074228 Environmental Express® K4300	Thermo Scientific™ Dionex™ ICS-2100, ICS-4000 Capillary HPIC, ICS-5000+ HPIC systems	Thermo Scientific™ Dionex™ AS-50 and AS-AP Autosampler	Thermo Scientific™ Dionex™ : 062374 (AS) 074938 (AS-AP) 069877 (AS-AP) 055056 (AS-50)