

April 13, 2018

RDE Sample Holder



MP-116 Sample holder is made to fit Spectroil M Series Rotating Disc Electrode Optical Emission Spectroscopy (RDEOES). Equivalent to Spectro Scientific M90200 & P-10524, SPC Science 070-010-101 and NSN 6640-01-445-1330. Made with Polypropylene copolymer and packaged 1000 cups/pack.

- Used for sample presentation to Spectroil M/N-W M/R instruments
- Can be used manually or with Double Disc Rapid Robot (D2R2)
- Used for ASTM 6595 and ASTM 6728 test methods